

# 3 FT ALUMINUM BRIEFCASE FOLDING RAMP

---

MR3B

321301



**RTITAN**  
RAMPS


## Operator's Manual



Read the Operator's Manual entirely. When you see this symbol, the subsequent instructions and warnings are serious follow without exception. Your life and the lives of others depend on it!

## PARTS DIAGRAM / EXPLODED VIEW

---



| KEY | BOX | DESCRIPTION           | QTY |
|-----|-----|-----------------------|-----|
| (1) | 1   | SECTION A             | 1   |
| (2) | 1   | SECTION B             | 1   |
| (3) | 1   | HINGE                 | 2   |
| (4) | 1   | CHAIN WITH FIXING PIN | 2   |
| (5) | 1   | HANDEL                | 1   |

## ASSEMBLY INSTRUCTIONS

---

- No assembly Required.

## ADDITIONAL INFORMATION

---

### SAFETY WARNING:

The maximum weight capacity is 500 LB. Do not exceed the maximum weight capacity.

The maximum loading height list is as follows:

#### **MAXIMUM \*OCCUPIED/\*\*UNOCCUPIED HEIGHT LIMIT TO GROUND LEVEL:**

2 FT Ramp: Occupied loading: 4-in maximum, unoccupied loading height: 6-in.

3 FT Ramp: Occupied loading: 6-in maximum, unoccupied loading height: 9-in.

4 FT Ramp: Occupied loading: 8-in maximum, unoccupied loading height: 12-in.

5 FT Ramp: Occupied loading: 10-in maximum, unoccupied loading height: 15-in.

6 FT Ramp: Occupied loading: 12-in maximum, unoccupied loading height: 18-in.

7 FT Ramp: Occupied loading: 14-in maximum, unoccupied loading height: 21-in.

8 FT Ramp: Occupied loading: 16-in maximum, unoccupied loading height: 24-in.

9 FT Ramp: Occupied loading: 18-in maximum, unoccupied loading height: 27-in.

10 FT Ramp: Occupied loading: 20-in maximum, unoccupied loading height: 30-in.

\*Occupied is a person sitting in the wheelchair/scooter while loading (2:12 slope or 9.6-degree angle)

\*\*Unoccupied is loading with nobody sitting in the wheelchair/scooter (3:12 slope or 14.5-degree angle)

#### **MAXIMUM WEIGHT CAPACITY IS 500 LB**

**(WEIGHT RATING IS BASED ON EVENLY DISTRIBUTED LOAD ON 2 AXLES)**

## **ACKNOWLEDGEMENT OF RISK AND RELEASE OF LIABILITY**

---

The use of any equipment, including this one, involves the potential risk of injury. Apart from any warranty claim that might be presented for a claimed defect in material or workmanship of the product, you accept and assume full responsibility for any and all injuries, damages (both economic and non-economic), and losses of any type, which may occur, and you fully and forever release and discharge Titan, its insurers, employees, officers, directors, associates, and agents from any and all claims, demands, damages, rights of action, or causes of action, present or future, whether the same be known or unknown, anticipated, or unanticipated, resulting from or arising out of the use of said equipment.

This equipment must be used with care by capable and competent individuals under supervision, if necessary.

## **TITAN LIMITED WARRANTY: TERMS, EXCLUSIONS AND LIMITATIONS OF REMEDIES**

---

This product comes with a one (1) year limited warranty that can be found at [www.titanramps.com/warranty.html](http://www.titanramps.com/warranty.html). Please review the same for all details regarding the Titan Limited Warranty.

**THE TITAN LIMITED WARRANTY FOUND AT [WWW.TITANRAMPS.COM/WARRANTY.HTML](http://WWW.TITANRAMPS.COM/WARRANTY.HTML) IS EXCLUSIVE AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY AND/OR FITNESS FOR A PARTICULAR PURPOSE, EACH OF WHICH IS HEREBY DISCLAIMED.**



**NEED HELP? CONTACT US FIRST.**

**1-800-605-7555**

[info@titanramps.com](mailto:info@titanramps.com)

[www.titanramps.com](http://www.titanramps.com)

© 2021 Titan Brands