

ALUMINUM THRESHOLD RAMP

TR2, TR3, TR4, TR5, TR6

323902, 323903, 323904, 323905, 323906



TITAN
RAMPS



Operator's Manual

Read the Operator's Manual entirely. When you see this symbol, the subsequent instructions and warnings are serious follow without exception. Your life and the lives of others depend on it!



KEY	DESCRIPTION	QTY
(1)	1" MODULAR THRESHOLD RAMP	1
(2)	2" MODULAR THRESHOLD RAMP EXT	1
(3)	MODULAR THRESHOLD RAMP PLATES	1
(4)	Screws FOR MODULAR THRESHOLD RAMP	3

MAXIMUM WEIGHT CAPACITY IS 600 LB

(WEIGHT RATING IS BASED ON EVENLY DISTRIBUTED LOAD ON 2 AXLES)

The ramp has a weight capacity of 300 lbs. with 1 axle and 600 lbs. with 2 axles

TR3 PARTS DIAGRAM / EXPLODED VIEW



KEY	DESCRIPTION	QTY
(1)	1" MODULAR THRESHOLD RAMP	1
(2)	2" MODULAR THRESHOLD RAMP EXT	1
(3)	3" MODULAR THRESHOLD RAMP EXT	1
(4)	MODULAR THRESHOLD RAMP PLATES	2
(5)	Screws FOR MODULAR THRESHOLD RAMP	5

**MAXIMUM WEIGHT CAPACITY IS 600 LB
(WEIGHT RATING IS BASED ON EVENLY DISTRIBUTED LOAD ON 2 AXLES)**

The ramp has a weight capacity of 300 lbs. with 1 axle and 600 lbs. with 2 axles



KEY	DESCRIPTION	QTY
(1)	1" MODULAR THRESHOLD RAMP	1
(2)	2" MODULAR THRESHOLD RAMP EXT	1
(3)	3" MODULAR THRESHOLD RAMP EXT	1
(4)	4" MODULAR THRESHOLD RAMP EXT	1
(5)	MODULAR THRESHOLD RAMP PLATES	3
(6)	SCREWS FOR MODULAR THRESHOLD RAMPS	7

MAXIMUM WEIGHT CAPACITY IS 600 LB

(WEIGHT RATING IS BASED ON EVENLY DISTRIBUTED LOAD ON 2 AXLES)

The ramp has a weight capacity of 300 lbs. with 1 axle and 600 lbs. with 2 axles

TR5 PARTS DIAGRAM / EXPLODED VIEW



KEY	DESCRIPTION	QTY
(1)	1" MODULAR THRESHOLD RAMP	1
(2)	2" MODULAR THRESHOLD RAMP EXT	1
(3)	3" MODULAR THRESHOLD RAMP EXT	1
(4)	4" MODULAR THRESHOLD RAMP EXT	1
(5)	5" MODULAR THRESHOLD RAMP EXT	1
(6)	MODULAR THRESHOLD RAMP PLATES	4
(7)	SCREWS FOR MODULAR THRESHOLD RAMPS	9

**MAXIMUM WEIGHT CAPACITY IS 600 LB
(WEIGHT RATING IS BASED ON EVENLY DISTRIBUTED LOAD ON 2 AXLES)**

The ramp has a weight capacity of 300 lbs. with 1 axle and 600 lbs. with 2 axles

TR6 PARTS DIAGRAM / EXPLODED VIEW



KEY	DESCRIPTION	QTY
(1)	1" MODULAR THRESHOLD RAMP	1
(2)	2" MODULAR THRESHOLD RAMP EXT	1
(3)	3" MODULAR THRESHOLD RAMP EXT	1
(4)	4" MODULAR THRESHOLD RAMP EXT	1
(5)	5" MODULAR THRESHOLD RAMP EXT	1
(6)	6" MODULAR THRESHOLD RAMP EXT	1
(7)	MODULAR THRESHOLD RAMP PLATES	5
(8)	SCREWS FOR MODULAR THRESHOLD RAMPS	11

**MAXIMUM WEIGHT CAPACITY IS 600 LB
(WEIGHT RATING IS BASED ON EVENLY DISTRIBUTED LOAD ON 2 AXLES)**

The ramp has a weight capacity of 300 lbs. with 1 axle and 600 lbs. with 2 axles

TR2 PARTS DIAGRAM / EXPLODED VIEW

MAXIMUM *OCCUPIED/UNOCCUPIED HEIGHT LIMIT TO GROUND LEVEL:**

- 2 FT Ramp: Occupied loading: 4-in maximum, unoccupied loading height: 6-in.
- 3 FT Ramp: Occupied loading: 6-in maximum, unoccupied loading height: 9-in.
- 4 FT Ramp: Occupied loading: 8-in maximum, unoccupied loading height: 12-in.
- 5 FT Ramp: Occupied loading: 10-in maximum, unoccupied loading height: 15-in.
- 6 FT Ramp: Occupied loading: 12-in maximum, unoccupied loading height: 18-in.
- 7 FT Ramp: Occupied loading: 14-in maximum, unoccupied loading height: 21-in.
- 8 FT Ramp: Occupied loading: 16-in maximum, unoccupied loading height: 24-in.
- 9 FT Ramp: Occupied loading: 18-in maximum, unoccupied loading height: 27-in.
- 10 FT Ramp: Occupied loading: 20-in maximum, unoccupied loading height: 30-in.

*Occupied is a person sitting in the wheelchair/scooter while loading (2:12 slope or 9.6-degree angle)

**Unoccupied is loading with nobody sitting in the wheelchair/scooter (3:12 slope or 14.5-degree angle)

MAXIMUM WEIGHT CAPACITY IS 600 LB

(WEIGHT RATING IS BASED ON EVENLY DISTRIBUTED LOAD ON 2 AXLES)

ACKNOWLEDGEMENT OF RISK AND RELEASE OF LIABILITY

The use of any equipment, including this one, involves the potential risk of injury. Apart from any warranty claim that might be presented for a claimed defect in material or workmanship of the product, you accept and assume full responsibility for any and all injuries, damages (both economic and non-economic), and losses of any type, which may occur, and you fully and forever release and discharge Titan, its insurers, employees, officers, directors, associates, and agents from any and all claims, demands, damages, rights of action, or causes of action, present or future, whether the same be known or unknown, anticipated, or unanticipated, resulting from or arising out of the use of said equipment.

This equipment must be used with care by capable and competent individuals under supervision, if necessary.

TITAN LIMITED WARRANTY: TERMS, EXCLUSIONS AND LIMITATIONS OF REMEDIES

This product comes with a one (1) year limited warranty that can be found at www.titanramps.com/warranty.html. Please review the same for all details regarding the Titan Limited Warranty.

THE TITAN LIMITED WARRANTY FOUND AT WWW.TITANRAMPS.COM/WARRANTY.HTML IS EXCLUSIVE AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY AND/OR FITNESS FOR A PARTICULAR PURPOSE, EACH OF WHICH IS HEREBY DISCLAIMED.



NEED HELP? CONTACT US FIRST.

1-800-605-7555

info@titanramps.com

www.titanramps.com

© 2021 Titan Brands